

Codonics Horizon Operating Software

Technical Brief

Version 4.3.1 Release Notes

Summary

Horizon® 4.3.1 software provides support for the Aeon motherboard and makes improvements to remote imager status handling.

Products Affected

GS/G2/G1 Imagers containing newer hardware identified by the serial number prefix 150C are affected. Please contact your Codonics representative for further information.

The Imager Software Version can be determined using the Control Panel by viewing the Imager Info screen from the Utilities menu.

NOTE: Horizon Imagers with serial number prefix 150C cannot be downgraded to any previous software version below 4.0.0. Imagers with a serial number prefix other than 150C cannot have their software upgraded to 4.0.0.

Product Improvements (Perfective Updates)

- ◆ Support for the Aeon Motherboard

- ◆ **Improved Error Message Logging.** Previous software releases logged more messages than necessary during remote status checks of Horizon via telnet. When checking status is done frequently this causes the logs to fill up quickly with unnecessary information. This effectively reduces the total amount of logging that can be maintained on a Horizon. Certain unnecessary messages are now no longer being logged during a remote status check of the imager.

Defects Corrected (Corrective Updates)

- ◆ **Updated the Linux software kernel.** When repeatedly checking status of the Horizon for use in a Kiosk via the Horizon SDK, occasional imager reboots can occur. Horizon 4.3.1 software uses an updated Linux kernel that corrects this issue.

Known Common Issues

- ◆ Previous issues can be found in previous Release Notes.

Known Uncommon Issues

- ◆ Previous issues can be found in previous Release Notes.

Technical Support

If problems occur during software installation, please contact Codonics Technical Support at any time.

Phone: +1.440.243.1198
Email: support@codonics.com
Website: www.codonics.com

Get it all with just one call
800.444.1198

All registered and unregistered trademarks are the property of their respective owners. Specifications subject to change without notice.
Patents: www.codonics.com/ip/patents
Copyright © 2018 Codonics, Inc. Printed in the U.S.A.
Part No. 901-462-009.01.



17991 Englewood Drive
Middleburg Heights, OH
44130 USA
+1.440.243.1198
+1.440.243.1334 Fax
Email info@codonics.com
www.codonics.com

Codonics Trading Co, Ltd.
317 Xianxia Rd. Building B
Unit 1412
Changning Dist., Shanghai
P.R. China, 200051
86-21-62787701
86-21-62787719 Fax

Codonics Limited KK
AQUACITY 9F,
4-16-23, Shibaura
Minato-ku, Tokyo,
108-0023 JAPAN
81-3-5730-2297
81-3-5730-2295 Fax