



Standard Model

Compact Model

## Intelligent Anesthesia Medication Station

Applicable Departments: Operating Room

Medicine Type: Toxic, Anesthetic, Psychoactive and Regular Drugs used during surgery





**Weight Detection**  
Gram Level Counting Algorithm



**<0.5s**  
Seconds Detection



**Steel Construction**  
SGCC Steel Plate, Highly Durable



**Intelligent Display**  
Individual Color Screen for each Bin with Smart Guidance



**Custom Configuration**  
Custom Layout as Needed



**Intelligent Alarm**  
Power Off/Forced Entry Alarm



**Voice Prompt**  
Automatic Error Warning



**Intelligent System**  
Dispense/Re-Supply/Inventory Check/Batch and Expiration Date Management



**Dual Authentication**  
Safety management of anesthetic medicine



**Dual Display**  
Reduce Errors, Ease of Operation



**Intelligent Lock Control**  
Authorized User Login, Multilevel Control

## Technical Specification

Physical Dimension	W*D*H: 700mm*695mm*1542mm W*D*H: 615mm*646mm*1491mm (Cabinet Only)
Material	Whole cabinet is made of high-strength aluminum alloy and metal with anodized surface
Intelligent Lock Control	Authorized user set up for cabinet lock access, allows individual set up for each bins
Intelligent Drug Bin Management	General, psychoactive, and anesthetic medicine control
Speaker	Voice Prompt
Display	Ten-point touch capacitive screen
Control Panel	Equipped with fingerprint reader, RFID reader, 7-inch touch screen display
Intelligent Monitoring	Intelligent real-time quantity and specification monitoring, accommodates oral, ampoule, vial, boxed, blister pack medicines
Intelligent Display	Each bin is equipped with a color LED display displaying the drug identification, name, specifications, quantity and operating instructions
Expandable Work Surface	Hidden, retractable
Casters	4 directional casters
Power Interface	AC, 100-240VAC, 300W
Network Port	Wired ethernet port
UPS	Backup UPS, 1000VA
Anti-theft Detection	Automatic alarm
Emergency Access	Dual mechanical locks ensure that the medicine can still be access when the power is down

### Optional Accessories

Empty bottle recycling bin, HD camera, prescription printer, barcode printer, storage box, catheter box, sharp object box

### Drawer Layout Top View

Standard Model				
Compact Model				