



CODONICS

Horizon[®] Imager

IHE Integration Statement

Part Number: 900-229-005, Rev A

Document / Software Version: 1.9.0

Date: 25 April 2006 [11:51]

This page intentionally left blank.

Revision	Date	Author	Description of Changes
1.9.0	24 Apr 2006	Jeremy Audino	Initial Release

IHE Integration Statement			
Vendor	Product Name	Version	Date
Codonics Inc.	Horizon Medical Imager (all models)	V1.9.0 and later	20 Jan 2006
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Consistent Presentation of Images	Print Server	N/A	
Internet address for vendor's IHE information: www.codonics.com/ihe			
Links to Standards Conformance Statements for the Implementation			
DICOM	www.codonics.com/Support/Horizon/dicom.php		
Links to general information on IHE			
In North America: www.rsna.org/IHE	In Europe: www.ihe-europe.org		In Japan: www.jira-net.or.jp/ihe-j