

# Codonics Color Paper / White Film

With more of today's medical systems using color applications, Codonics leads the industry with our ChromaVista Color Paper





"We work with a neurology group that was thrilled when they saw how our razor-sharp color studies on paper could be slipped conveniently into patient charts. It literally changed their practice."

*John C. Litchney, MBA, RT, Administrator for The Cleveland Clinic, Department of Regional Radiology*

## A4 Size

## A Size



- Codonics exclusive Medical Color Matching (MCM™) ensures prints match soft-copy monitors
- Ideal for referring physicians and patient take-homes
- Specifically designed for medical applications, Codonics color paper is FDA approved and fully archivable
- Image quality is far superior to laser and office printers featuring 16.7 million colors
- The only color output in the industry that truly rivals conventional color photographic film



**CODONICS®**


We bring the future into focus

# Codonics ChromaVista® Color Paper / White Film



Codonics introduced the world's first DICOM color imager and has been pioneering multi-media imagers for over two decades. We are the industry leader in color imagers with nearly 30,000 installations worldwide.

Our Horizon Multi-media Imager features stunning edge-to-edge color printing, maximizing imaging size and reducing waste. Codonics exclusive Medical Color Matching (MCM™) is an advanced image processing feature that adjusts printed output colors to accurately match soft-copy monitors. This process corrects differences in hue, saturation, and intensity between the imager and soft-copy review stations. User's can easily determine the best match for their screen using the automatic MCM bracketing feature which is built into every Codonics imager.



Color prints always match softcopy monitors with our exclusive Medical Color Matching (MCM™) feature. Consistency, clarity and saturation with Codonics color dye-diffusion technology far exceeds that of wax or inkjet printers.


ChromaVista Color Paper / White Film has a heavy, photo-weight base material that resists water and tearing. Since there is no ink or toner, images never smear or transfer onto hands or other prints.

## Specifications:

- Compatible with the Horizon family of multi-media imagers
- Continuous tone color rendering
- $D_{min} < 0.10$  OD,  $D_{max} > 2.50$  OD (reflective)
- Thickness ~ 0.2mm

## Available Sizes:

- A (8.5" x 11")
- A4 (210mm x 297mm)



**CODONICS®**

We bring the future into focus

Call Codonics today at **800.444.1198** or visit [www.codonics.com](http://www.codonics.com) for more information.

All registered and unregistered trademarks are the property of their respective owners. Specifications subject to change without notice. Patents: [www.codonics.com/ip/patents](http://www.codonics.com/ip/patents) Copyright©2006-2022 Codonics Inc. 02/2022



17991 Englewood Drive  
Middleburg Heights, OH 44130 USA  
+1.440.243.1198  
+1.440.243.1334 Fax  
Email [info@codonics.com](mailto:info@codonics.com)  
[www.codonics.com](http://www.codonics.com)